





Local Structure Analysis of Non-Crystalline Materials based on PDFs

Article by Akito Sasaki, Rigaku Corporation
Presented at Denver X-ray Conference 2019

Experimental and analysis details


Using the SmartLab system, we obtained high-quality scattering patterns for PDF analysis. We used a high-power PhotonMax 9 kW rotating anode X-ray source with Ag target and CBO-E elliptical multilayer mirror for Ag radiation to condition a high-resolution focusing beam. A scattering protector and a 2.5-deg Soller slit were used in the receiving optics to eliminate parasitic scattering from the incident optics and air. We also used a D/teX Ultra HE 1D detector as a high-efficiency detector for the high-energy X-ray source to collect data rapidly.



SmartLab Studio II PDF plugin


Experimental results in comparison with synchrotron results

We used the Ag K α wavelength source to collect as high-q data as possible. Compared to the structure factors from the experimental results from a synchrotron facility, we obtained almost the same $S(Q)$ pattern using SmartLab. As a result, the calculated PDF based on the results from SmartLab were quite similar to that from a synchrotron facility.



Structure analysis of SiO₂ and its derivatives based on experimental PDFs


Using the SmartLab system, we collected scattering patterns of noncrystalline SiO₂ and its derivatives (77SiO₂+23Na₂O, 77SiO₂+23K₂O). Based on the scattering patterns, the PDF of each compound was calculated. Based on the obtained PDFs, the structure of each compound was analyzed using a reverse Monte Carlo algorithm (RMC Profile software¹). Each structural model well reproduced the corresponding experimentally obtained PDF.



Local structure of rings in crystalline (left) and amorphous silica (right)

Then, the numbers of n -membered rings existing in each structure were counted using the ISAACS software^{2,3}.

In the case of SiO₂, 12-membered rings were the most probable⁴ and rings larger than 20-membered were not observed. In the cases of Na,K+SiO₂, 8-membered rings were the most probable and rings larger than 20-membered were also observed for both cases. This indicates that Si/O rings became larger due to the incorporation of sodium/potassium ions into them.



References

- (1) R. L. McGreevy *et al.*, *Molecular Simulation*, **1**, 359-367(1988).
- (2) S. Le Roux and V. Petkov: *J. Appl. Cryst.*, **43** (2010), 181-185.
- (3) L. Guttman: *J. Non-Cryst. Solids*, **116** (1990), 145-147.
- (4) S.Kohara *et al.*, *Proc. Nat.Acad.Sci.*, **108(36)**, 14780(2011).